


Prepare Scan

1. Fill in the necessary user fields. For the **date of birth**, use a **specific format** of **MM-DD-YYYY** separated by hyphens.

e.g. 01-30-2005 03-01-2005

2. Click



. This will automatically clear the fields

Test Placement and Signals

3. Click



on the **upper-left corner** to test the real-time plot

Start Recording (64 Seconds)

4. Click



if there exists an existing test to reset the test

5. **Simultaneously** hold the



key down on the keyboard, navigate the mouse

cursor above



on the upper-left corner and click the button.

Scan Finishes

Time Span Left: 000d, 00:00:58.800



Scan'g

Logger on


PLY

SET


(Scan finishes if time runs down to 0)

Save to CSV File


6. On the upper-left corner:



Click: **File**




Click: **Save As ...**




You will be prompted with this window.

Click on the **Save as type:** row




Click the **CSV** option

7. In the file name field, type without blank spaces, just the number: 1




8. DO NOT CLICK:



Instead, click the




key.



Post-Scan Plot

9. In the window with the fields for patient information, click

Plot Data



10. Click

Quit

to exit the window.

BACK TO STEP 1 FOR NEW PATIENT

Quit Program (If Necessary)

11. Click



key, then close the patient field window.

Enter Fields

Initials: Date of Birth (MM-DD-YYYY): Sex: ☒ Male ☐ Female

Save Fields Plot Data